

# ALB Anti Serum

**Product Name:** Goat Anti Human Albumin, Serum

**Product Cat#:** BBMALBN102

**Source:** Goat

**Antibody Type:** Anti Serum

**Assay:**  $\geq 50\text{mg/ml}$

**valence of antibody:** 1: 5

**Specificity:** Humanized ALB

**Buffer:** 20mM Tris+250mM NaCl+15mM NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Storage:** 2~8°C

**Validity:** 3 years

**Application:** Immunoturbidimetry

**Reagent Preparation:**

Reagent 2 buffer: 20mM Tris+250mM NaCl+15mM NaN<sub>3</sub>, pH7.4

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 3 ~ 1: 10

*The above is the recommended formula and usage, customers can adjust according to their needs*