

# ApoA I Anti Serum

**Product Name:** Goat Anti Human Apolipoprotein A<sub>2</sub>, Serum

**Product Cat#:** BBNAPOA1N101

**Source:** Goat

**Antibody Type:** Anti Serum

**Assay:**  $\geq 30$ mg/ml

**valence of antibody:** 1:5

**Specificity:** Humanized ApoA I

**Buffer:** 20mM Tris +250mM NaCl+0.09%NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Application:** Immunoturbidimetry

**Reagent Preparation:**

Reagent 1: 50mM Tris+250mM NaCl+0.1%EDTA2Na+0.1%NaN<sub>3</sub>+5%PEG6000+0.5%Tween20, pH7.4

Reagent 2 buffer: 50mM Tris+100mM NaCl+0.1%EDTA2Na+0.1%NaN<sub>3</sub> +0.2%Tween20, pH7.4

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 5

*The above is the recommended formula and usage, customers can adjust according to their needs*