

ASO Positive Serum Calibrator, High Level

Product Name: ASO Positive Serum Calibrator, High Level

Product Cat#: FPZ0174-1

Source: Goat

Molecular weight: 150kDa

Assay: ≥ 5000 IU/ml

Specificity: Humanized ASO

Buffer: 20mM Tris + 150mM NaCl + 0.09% NaN₃, pH7.4

Format: Liquid

Storage: -20 ± 5 °C

Validity: 3 years

Application: Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

The above is the recommended formula and usage, customers can adjust according to their needs