

## Monoclonal Mouse Anti-Human CRP

**Product Name:** Monoclonal Mouse Anti-Human CRP, Affinity Purified

**Product Cat#:** BBNCRPN102

**Source:** Mouse

**Antibody Type:** Affinity Purified

**Assay:**  $\geq 1\text{mg/ml}$

**Purity:**  $> 90\%$

**Subtype:** IgG1

**Molecular mass:** 150kDa

**Specificity:** Anti-Human CRP

**Buffer:** 10mM PB + 150mM NaCl + 0.1% P300, pH7.4

**Format:** Liquid

**Storage:**  $-20 \pm 5^\circ\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

**Attention:**

1. The amount of subpackage can be selected at will.
2. This product is a biologically active protein, which is only used for research and development and production purposes, and cannot be directly used for diagnosis purposes.
3. Pay attention to storage conditions.