

## Apo AI Calibrator, High Level

**Product Name:** Apo AI Calibrator, High Level

**Product Cat#:** GBNAPOA101

**Source:** Recombinant proteins

**Molecular weight:** 28kDa

**Assay:**  $\geq 2$ mg/ml

**Specificity:** Humanized Apo AI

**Buffer:** 50mM Tris + 150mM NaCl + 0.05% NaN<sub>3</sub>, pH8.0

**Format:** Liquid

**Storage:**  $-20 \pm 5^{\circ}\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*