

Apo B Calibrator, High Level

Product Name: Apo B Calibrator, High Level

Product Cat#: GBNAPOB101

Source: Recombinant proteins

Molecular weight: 550kDa

Assay: $\geq 2\text{mg/ml}$

Specificity: Humanized Apo B

Buffer: 50mM Tris + 150mM NaCl + 0.05% NaN₃, pH8.5

Format: Liquid

Storage: $-20 \pm 5^\circ\text{C}$

Validity: 3 years

Application: Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

The above is the recommended formula and usage, customers can adjust according to their needs