

## Certificate of Analysis

**Product name:** Anti-MYO polyclonal antibody

**Catalog number:** BBMYON101

**Expression system:** Rabbit

**Appearance:** Liquid

**Purity:** /

**Isotype:** /

**Molecular weight:** 150 kDa

**Clone:** /

**Storage buffer:** 10mM PB+100mM NaCl+0.09% NaN<sub>3</sub>, pH7.4

**Storage condition:** -20±5°C

**Shelf life:** 3 years

**Intended use:** /

**Application:** Immunoturbidimetry

### Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.