

Certificate of Analysis

Product name: ApoB quality control (Goat serum)

Catalog number: BBNAPON102

Expression system: Goat

Appearance: Liquid

Purity: /

Isotype: /

Molecular weight: 150 kDa

Clone: /

Storage buffer: 25mM Tris+250mM NaCl+0.09% NaN₃, pH8.0

Storage condition: -20±5°C

Shelf life: 3 years

Intended use: /

Application: Immunoturbidimetry

Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.