

Datasheet

Product name: Goat Anti Human C-Reactive Protein Polyclonal Antibody

Catalog number: BBNCRPN104

Expression system: Goat

Appearance: Liquid

Purity: /

Isotype: /

Molecular weight: 150 kDa

Clone: /

Storage buffer: 10mM Tris+150mM NaCl+5mM EDTA+0.1% NaN₃, pH8.0

Storage condition: -20±5°C

Shelf life: 3 years

Intended use: Anti-Human C-Reactive Protein

Application: Immunoturbidimetry

Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.