

Datasheet

Product name: Anti-CKMB monoclonal antibody

Catalog number: BRJCKMBS101

Expression system: Mouse

Appearance: Liquid

Purity: $\geq 90\%$

Isotype: IgG1

Molecular weight: 150 kDa

Clone: /

Storage buffer: 10mM PB+150mM NaCl+0.1% P300, pH7.4

Storage condition: $-20\pm 5^{\circ}\text{C}$

Shelf life: 3 years

Intended use: Anti-Human Creatine Kinase-MB Isoenzyme Protein

Application: Immunoturbidimetry

Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.