

# CRP Calibrator, High Level

**Product Name:** CRP Calibrator, High Level

**Product Cat#:** GBNCRPN101

**Source:** Natural extraction

**Molecular weight:** 120kDa

**Assay:** 2mg/ml

**Specificity:** Humanized CRP

**Buffer:** 20mM Tris + 280mM NaCl + 5mM CaCl<sub>2</sub> + 0.1% NaN<sub>3</sub>, pH8.0

**Format:** Liquid

**Storage:** 2~8°C

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*