

CysC Calibrator, High Level

Product Name: CysC Calibrator, High Level

Product Cat#: GBNCYSN101

Source: Recombinant proteins

Molecular weight: 13kDa

Assay: ≥ 1.5 mg/ml

Specificity: Humanized CysC

Buffer: 50mM Tris + 150mM NaCl + 3mM EDTA, pH8.5

Format: Liquid

Storage: -20 ± 5 °C

Validity: 3 years

Application: Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

The above is the recommended formula and usage, customers can adjust according to their needs