

# D-dimer Calibrator, High Level

**Product Name:** D-dimer Calibrator, High Level

**Product Cat#:** GBNDIMN101

**Source:** Recombinant proteins

**Molecular weight:** 190kDa

**Assay:**  $\geq 1$ mg/ml

**Specificity:** Humanized D-dimer

**Buffer:** 10mM PBS + 150mM NaCl + 0.09% NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Storage:**  $-20 \pm 5^{\circ}\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*