

# Datasheet

**Product name:** RF Antigen

**Catalog number:** GBNIGGN101

**Expression system:** /

**Appearance:** Liquid

**Purity:** /

**Molecular weight:** 150 kDa

**Storage buffer:** 10mM PB +150mM NaCl +0.1% NaN<sub>3</sub>, pH7.4

**Tag:** /

**Storage condition:** -20±5°C

**Shelf life:** 3 years

**Intended use:** /

**Application:** TIA/LETIA

## Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.