

Datasheet

Product name: Streptolysin O (SLO) Recombinant Antigen

Catalog number: GBNSLON106

Expression system: Recombinant protein

Appearance: Liquid

Purity: /

Molecular weight: 450 kDa

Storage buffer: 20mM PB + 150mM NaCl + 0.05% NaN₃, pH7.4

Tag: /

Storage condition: -20±5°C

Shelf life: 3 years

Intended use: /

Application: LETIA

Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.