

IgA/C3/C4 Calibrator, High Level

Product Name: IgA/C3/C4 Calibrator, High Level

Product Cat#: FPZ0728

Source: Recombinant proteins

Molecular weight: /

Assay: IgG \geq 45g/L; IgA \geq 8g/L; IgM \geq 4g/L; C3 \geq 4g/L; C4 \geq 1g/L; TRF \geq 4g/L

Specificity: Humanized

Buffer: 10mM PBS + 150mM NaCl + 0.09% NaN₃, pH7.4

Format: Liquid

Storage: -20 \pm 5 °C

Validity: 3 years

Application: Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

The above is the recommended formula and usage, customers can adjust according to their needs