

# IgA Anti Serum

**Product Name:** Goat Anti Human IgA ( $\alpha$  Chain Specific), Serum

**Product Cat#:** BBNIGAN102

**Source:** Goat

**Antibody Type:** Anti Serum

**Assay:**  $\geq 30\text{mg/ml}$

**valence of antibody:** 1: 4

**Specificity:** Humanized IgA

**Buffer:** 25mM Tris + 250mM NaCl, pH7.5

**Format:** Liquid

**Storage:**  $-20 \pm 5^\circ\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

**Reagent Preparation:**

Reagent 2 buffer: 25mM Tris + 250mM NaCl, pH7.5

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 3 ~ 1: 10

*The above is the recommended formula and usage, customers can adjust according to their needs*