

IgG Anti Serum

Product Name: Goat Anti Human IgG (high activity), Serum

Product Cat#: BBNIGGN104

Source: Goat

Antibody Type: Anti Serum

Assay: $\geq 30\text{mg/ml}$

valence of antibody: 1: 4

Specificity: Humanized IgG

Buffer: 20mM Tris+250mM NaCl+0.09% NaN₃, pH7.4

Format: Liquid

Storage: $-20 \pm 5^{\circ}\text{C}$

Validity: 3 years

Application: Immunoturbidimetry

Reagent Preparation:

Reagent 2 buffer: 20mM Tris+250mM NaCl+0.09% NaN₃, pH7.4

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 3 ~ 1: 10

The above is the recommended formula and usage, customers can adjust according to their needs