

# IgM Calibrator, High Level

**Product Name:** IgM Calibrator, High Level

**Product Cat#:** GBNIGMN102

**Source:** Recombinant proteins

**Molecular weight:** 970kDa

**Assay:**  $\geq 4\text{g/L}$

**Specificity:** Humanized IgM

**Buffer:** 10mM PBS + 150mM NaCl + 0.09% NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Storage:**  $-20 \pm 5^\circ\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*