

IgM Anti Serum

Product Name: Goat Anti Human IgM (μ Chain Specific),Serum

Product Cat#: BBNIGMN102

Source: Goat

Antibody Type: Anti Serum

Assay: $\geq 30\text{mg/ml}$

valence of antibody: 1: 4

Specificity: Humanized IgM

Buffer: 50mM Tris + 500mM NaCl + 0.09%NaN₃, pH7.5

Format: Liquid

Storage: $-20 \pm 5^\circ\text{C}$

Validity: 3 years

Application: Immunoturbidimetry

Reagent Preparation:

Reagent 2 buffer: 50mM Tris + 500mM NaCl + 0.09%NaN₃, pH7.5

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 3 ~ 1: 10

The above is the recommended formula and usage, customers can adjust according to their needs