

## Lp(a) Calibrator, High Level

**Product Name:** Lp(a) Calibrator, High Level

**Product Cat#:** GBNLPA103

**Source:** Recombinant proteins

**Molecular weight:** 55kDa

**Assay:**  $\geq 1\text{mg/ml}$

**Specificity:** Humanized Lp(a)

**Buffer:** 20mM Tris + 200mM NaCl + 0.05% NaN<sub>3</sub>, pH8.0

**Format:** Liquid

**Storage:**  $-20 \pm 5^\circ\text{C}$

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*