

## RF Calibrator, High Level

**Product Name:** RF Calibrator, High Level

**Product Cat#:** FPZ0432

**Source:** Recombinant expression

**Molecular weight:** 150kDa

**Assay:** Based on specific batches

**Specificity:** Humanized RF

**Buffer:** 20mM Tris-HCl +0.09% NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Storage:** -20±5℃

**Validity:** 3 years

**Application:** Immunoturbidimetry

Using customer's self-matching buffer or the matching buffer provided by the manufacturer diluted to the appropriate concentration and then used

*The above is the recommended formula and usage, customers can adjust according to their needs*