

Rabbit Anti Human C-Reactive Protein Antibody

Product Name: Rabbit Anti Human C-Reactive Protein, Affinity Purified Antibody

Product Cat#: BBNCRPN101

Source: Rabbit

Antibody Type: Affinity Purified Antibody

Assay: $\geq 1\text{mg/ml}$

valence of antibody: /

Specificity: Humanized CRP

Buffer: 10mM PBS + 150mM NaCl + 0.1% P300, pH8.0

Format: Liquid

Storage: 2~8°C

Validity: 3 years

Application: Immunoturbidimetry

Attention:

1. The amount of subpackage can be selected at will.
2. This product is a biologically active protein, which is only used for research and development and production purposes, and cannot be directly used for diagnosis purposes.
3. Pay attention to storage conditions.