

# ApoB Anti Serum

**Product Name:** Goat Anti Human Apolipoprotein B, Serum

**Product Cat#:** BBNAPON101

**Source:** Goat

**Antibody Type:** Anti Serum

**Assay:**  $\geq 30\text{mg/ml}$

**valence of antibody:** 1: 6

**Specificity:** Humanized ApoB

**Buffer:** 25mM Tris +250mM NaCl+0.09%NaN<sub>3</sub>, pH7.4

**Format:** Liquid

**Application:** Immunoturbidimetry

**Reagent Preparation:**

Reagent 1: 50mM Tris+250mM NaCl+0.1%EDTA2Na+0.1%NaN<sub>3</sub>+5%PEG6000+0.5%Tween20, pH7.4

Reagent 2 buffer: 50mM Tris+100mM NaCl+0.1%EDTA2Na+0.1%NaN<sub>3</sub> +0.2%Tween20, pH7.4

Anti Serum should be diluted by reagent 2 buffer at ratio 1: 6

*The above is the recommended formula and usage, customers can adjust according to their needs*