

Datasheet

Product name: Rabbit Anti Human β -2-Microglobulin, Affinity Purified Antibody

Catalog number: BBNB2MN101

Expression system: Recombinant protein

Appearance: liquid

Purity: /

Isotype: /

Molecular weight: 150 kDa

Clone: /

Storage buffer: 10mM PB + 150mM NaCl + 0.05% NaN₃, pH7.4

Storage condition: -20 \pm 5 $^{\circ}$ C

Shelf life: 3 years

Intended use: /

Application: Immunoturbidimetry

Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.