

## Datasheet

**Product name:** Anti-CKMB monoclonal antibody

**Catalog number:** BRJCKMBS102

**Expression system:** Mouse

**Appearance:** Liquid

**Purity:**  $\geq 90\%$

**Isotype:** IgG1

**Molecular weight:** 150 kDa

**Clone:** /

**Storage buffer:** 10mM PB+150mM NaCl+0.1% P300, pH7.4

**Storage condition:**  $-20\pm 5^{\circ}\text{C}$

**Shelf life:** 3 years

**Intended use:** Anti-Human Creatine Kinase-MB Isoenzyme Protein

**Application:** Immunoturbidimetry

### Notes:

1. This product is a biologically active protein intended for use in research and production, not for use in diagnostic procedures.
2. Store this product under the specified conditions.