

Rabbit Anti Human Ferritin Protein

Product Name: Rabbit Anti Human Ferritin Protein, Affinity Purified Antibody

Product Cat#: BBNFERN101

Source: Rabbit

Antibody Type: Affinity Purified Antibody

Assay: $\geq 1\text{mg/ml}$

valence of antibody: /

Specificity: Humanized Ferritin Protein

Buffer: 20mM Heps+300mM NaCl+0.09% NaN₃, pH 7.5

Format: Liquid

Storage: $-20 \pm 5^{\circ}\text{C}$

Validity: 3 years

Application: Immunoturbidimetry

Reagent Preparation:

Latex enhanced immune turbidity platform, coated with latex.